








SETON RE-SEAL CROWNS

Ultra

707-782-0900
 707-782-0902 fax
sales@setonscientific.com
 P.O. Box 750003
 Petaluma, CA 94975 USA

ALUMINUM CROWNS

Rev.2

 <p>4065 +Aluminum, clear anodize +60 KRPM Max +Use Tubes: 9053 +For Rotors: 70Ti,60Ti,55.2Ti, 50.2Ti</p>	 <p>4068 +Aluminum, clear anodized +75 KRPM Max +Use Tubes: 9041R +For Rotors: 90Ti,80Ti,75Ti, 70.1Ti</p>	 <p>4069 +Aluminum, Black anodized +65 KRPM Max +Use Tubes: 9041R +For Rotors: Ty65,50Ti</p>
 <p>4491 +Aluminum, Black anodized +80 KRPM Max +Use Tubes: 9042, 9043 +For Rotors: S80-AT3,S58-A</p>	 <p>4257 +Aluminum, Black anodized +100 KRPM Max +Use Tubes: 9028 +For Rotors: S100-AT6</p>	 <p>4146 +Aluminum, clear anodize +60 KRPM Max +Use Tubes: 9054 +For Rotors: S55-A</p>
 <p>4492 +Aluminum, clear anodize +45 KRPM Max +Use Tubes: 9060 +For Rotors: 45Ti</p>		