

SETON

CONE-TOP RE-SEAL TUBES Polyallomer Thin Wall Ultra

707-782-0900
707-782-0902 fax
sales@setonscientific.com
P.O. Box 750003
Petaluma, CA 94975 USA

CONE-TOP RE-SEAL PA TUBES

Rev. 2

Tube Size	Seton PN Box of 50	Volume Nom/Actual ml	Use Crown (Plug)(1)	k-factor (Rotor)	For Beckman Rotor	For Sorvall (Thermo) Rotor	Like Beckman Tube PN(2)	Like Sorvall Tube PN
7/16x1" (11x25mm)	9011	1.5/1.5	4246 (4075)	14 (TLA100.2)	TLS55, TLA-120.2 MLA-130	See (5)	344624	75000501
7/16x1-1/4" (11x32mm)	9012	2.0/1.9	4246 (4075)	21 (SW60)	TLS55	S55-S, TH660, TFT80.4	344625	75000459
7/16x2-3/8" (11x60mm)	9059	3.4/3.4	4246 (4075)		SW60	TH660 TFT80.4	---	75000483
1/2x1" (13x25mm)	9024	2.0/1.9	4247 (4075)	18 (SW65)	SW55, 100Ti 50.4Ti	TFT45.6 AH-650	345829	75000458
1/2x1-1/4" (13x32mm)	9025	3.5/3.1	4247 (4075)	27 (50.3Ti)	SW55	TFT45.6 AH-650	349621	75000456
1/2x2" (13x48mm)	9058	4.0	4247 (4075)	12 (TLA100.4)	SW55	TFT45.6 AH-650	---	75000482
1/2x2" (13x51mm)	9029	5.0/5.1	4249 (4075)		TLA100.4	---	362248	---
9/16x1-7/8" (14x47mm)	9032	5.9/6.1	4219 (4075)	61 (SW41)	SW41 SW40	TH-641	355537	75000454
9/16x3-1/2" (14x88mm)	9055	10/10.2	4219 (4075)		SW41 SW40	TH-641	---	75000411
5/8x1-3/4" (16x45mm)	9040	6.3/6.8	4248 (4075)	98 ((SW28.1)	See (3)	See (6)	345830	75000450
5/8x4" (16x100mm)	9056	15.0	4248 (4075)		SW32.1 SW28.1	Surespin 630/17 AH-629/17	356291	75000414
1x1-1/2" (25x38mm)	9050	15/14.7	4216 (4075)	102 (SW28)	See (4)	Surespin 630/36 AH-629/36	343664	---
1x3-1/4" (25x81mm)	9051	30/29.7	4216 (4075)		SW32 SW28	Surespin 630/36 AH-629/36	344623	7500427

©2000-2015 Seton Scientific

NOTES:

- (1) PN 4075 is a package of ten plugs; a single plug is PN 4073.
- (2) Beckman Bell-Top and Seton Cone-Top tubes are interchangeable. Seton Cone-Top tubes can also be used with Beckman Bell-Top spacers.
- (3) Used in the following rotors: 90Ti, 80Ti, 75Ti, 70.1Ti, 65, 50Ti, 50, 40., SW32.1 and SW28.1.
- (4) Used in the following rotors: 70Ti, 60Ti, 55.2Ti, 50.2Ti, SW32 and SW28.
- (5) Used in the following rotors: TH660, TFT80.4., S55-S, S150-AT, S140-AT and S120-AT2.
- (6) Used in the following rotors: T-890, T-1270, Surespin 630/17 and AH-629/17.