







SETON

CONE-TOP CROWNS Re-Seal *Ultra*

707-782-0900
707-782-0902 fax
sales@setonscientific.com
P.O. Box 750003
Petaluma, CA 94975 USA

NORYL & DELRIN CROWNS

Rev.2

 <p>4216 +Noryl +70 KRPM Max +Use Tubes: 9050, 9051 +For Rotors: 70Ti,60Ti,55.2Ti, 50.2Ti,SW32,SW28</p>	 <p>4219 +Noryl +41 KRPM Max +Use Tubes: 9032, 9055 +For Rotors: 90Ti,80Ti,75Ti, 70.1Ti</p>	 <p>4246 +Noryl +100 KRPM Max +Use Tubes: 9011, 9012 +For Rotors: TLS55, S55-S, TH660,TFT80.4</p>
 <p>4247 +Noryl +100 KRPM Max +Use Tubes: 9024, 9025, 9058 +For Rotors: SW55,100Ti,50.4Ti TFT45.6,AH-650</p>	 <p>4248 +Noryl +80 KRPM Max +Use Tubes: 9040, 9056 +For Rotors: 90Ti,80Ti,75Ti,65 50Ti,50,40,SW32.1 SW28.1</p>	 <p>4249 +White Delrin +100 KRPM Max +Use Tubes: 9029 +For Rotors: TLA-100.4</p>

©2000-2015 Seton Scientific