

SETON

SELF-SEAL TUBES & CROWNS Thin Wall Polyallomer *Ultra*

707-782-0900
707-782-0902 fax
sales@setonscientific.com
P.O. Box 750003
Petaluma, CA 94975 USA

SELF-SEAL TUBES & CROWNS

Rev.2

Tube Size	Seton PN Box of 50	Tube Capacity ml	Cap Assy PN Seton (Thermo Sorvall)	For Thermo-Sorvall Rotor	For Kontron Rotor	Like Sorvall Tube PN	Max Speed KRPM
5/16x2" (8x51mm)	5009	1.5	4510 (54459)	TFT80.2	TFT80.2	54460	80
7/16x2-3/8" (11x60mm)	5012	4.4	4447 (54366)	TFT80.4	TFT80.4	03955	80

©2000-2015 Seton Scientific