













SETON EASY-SEAL CROWNS

Ultra

707-782-0900
 707-782-0902 fax
sales@setonscientific.com
 P.O. Box 750003
 Petaluma, CA 94975 USA

NORYL, DELRIN & ALUMINUM CROWNS

Rev.3

 <p>4015 Assy +Black Noryl +60 KRPM Max +Use Tubes: 6052, 8052 +For Rotors: 70Ti,60Ti,55.2Ti, 50.2Ti,42.1 +Hex Driver 4013</p>	 <p>4016 Assy +White Delrin +60 KRPM Max +Use Tubes: 6052, 8052 +For Rotors: 70Ti,60Ti,55.2Ti, 50.2Ti,42.1 +Hex Driver 4013</p>	 <p>4017 Assy +Black Noryl +65 KRPM Max +Use Tubes: 6041, 8041 +For Rotors: 90Ti,80Ti, 75Ti,70.1Ti +Hex Driver 4013</p>
 <p>4018 Assy +White Delrin +65 KRPM Max +Use Tubes: 6041, 8041 +For Rotors: 90Ti,80Ti, 75Ti,70.1Ti +Hex Driver 4013</p>	 <p>4019 Assy +Black Noryl +65 KRPM Max +Use Tubes: 6041, 8041 +For Rotors: 65,50Ti,40 A-1256,T-1270 +Hex Driver 4013</p>	 <p>4020 Assy +White Delrin +65 KRPM Max +Use Tubes: 6041, 8041 +For Rotors: 65,50Ti,40 A-1256,T-1270 +Hex Driver 4013</p>
 <p>4063 & 4064 Assy +Noryl, Delrin, Aluminum +80 KRPM Max +Use Tubes: 6041, 8041 +For Rotors: 90Ti,80Ti, 75Ti,70.1Ti +Hex Driver 4013 +Vise 4184 (1)</p>	 <p>4171 Assy +Delrin, Aluminum +50 KRPM Max +Use Tubes: 6023, 8023 +For Rotors: 50.4Ti,50.3Ti,40.2 40.3,TFT-45.6 +Hex Driver 4189</p>	 <p>4264 Assy +Delrin, Aluminum +120 KRPM Max +Use Tubes: 6011, 8011 +For Rotors: ?? ?? +Hex Driver 4189 +Vise 4172</p>
 <p>4267 Assy +Delrin, Aluminum +120 KRPM Max +Use Tubes: 6011, 8011 +For Rotors: TLA100.3, TLA100.4 S100-AT6, S100-AT4 +Hex Driver 4189 +Vise 4172</p>	 <p>4202 Assy +White Delrin +45 KRPM Max +Use Tubes: 6060, 8060 +For Rotors: 45Ti, Ty35 +Hex Driver 4013</p>	 <p>4111 Assy +Black Noryl +32 KRPM Max +Use Tubes: 6051 +For Rotors: SW32, SW28, SW +Hex Driver 4013</p>

©2000-2016 Seton Scientific

(1) Also can use Beckman vise PN 305075.