

707-782-0900
707-782-0902 fax
sales@setonscientific.com
P.O. Box 750003
Petaluma, CA 94975 USA

Page 1/2

OSR THREADED DELRIN TUBES

Rev. 2

Tube Size	Volume ml Nom/Fill	Assembly PN (1)	PN Tube (4)	PN Cap Assy (2)	PN O-ring (Viton)	Max Speed KRPM Fill/Part Fill(3) 4 deg C	Max Speed KRPM Fill/Part Fill(3) 25 deg C	For Beckman Rotors	For Sorvall/Kontron Rotors
7/16x2-3/8" (11x60mm)	2 2/21	4522	4532	4389	4390	1 1/10	1 1/9	---	TFT80.4
1/2x2-1/2" (13x64mm)	2 1/6	4523	4533	4279	4280	1	1 1/8	50.4Ti, 50.3Ti 40.3, 40.2	TFT45.6
5/8x3" (16x76mm)	1 1/2	4524	4534	4283	4284	1 1/4	1	see (5)	T-1270, T-865, A-1256, TFT80.13
1x3-1/2" (25x89mm)	1 27/50	4525	4535	4287	4288	1 3/7	1 3/7	S6e (5)	TFT: 70.38, 65.38, 55.38, 5 0.38
1-1/2x4" (38x102mm)	1 22/25	4516	4517	4291	4292	1 1/7	1 1/7	45Ti, 35	TFT45.94

©2000-2015 Seton Scientific

- (1) Assembly consists of one Delrin tube, aluminum cap and Viton® o-ring
- (2) Assembly consists of one aluminum cap and Viton o-ring.
- (3) Tubes must be centrifuged with caps and tubes must be full to use at any speed.
- (4) Delrin tubes must *not* be autoclaved.
- (5) Used in the following rotors: 90Ti, 80Ti, 75Ti, 70.1Ti, 65, 50Ti, 50 & 40.
- (6) Used in the following rotors: 70Ti, 60Ti, 55.2Ti, 50.2Ti, 42.1. Maximum speed is 40KRPM in the Beckman 50.2Ti rotor.

OSR TUBES FOR USE IN SW ROTORS

Rev 2

Tube Size	Volume ml Nom/Fill	PN Tube	Tube Material	Max Speed KRPM (3) 4 deg C	For Beckman Rotors	For Sorvall/Kontron Rotors
7/16x2-3/8" (11x60mm)	4.4/2.1	4519	Delrin(1)	45	SW60Ti	TH-660
		4528	Kynar(2)			
1/2x2" (13x51mm)	6.5/2.4	4520	Delrin	45	SW55Ti	AH-650
		4529	Kynar			
9/16x3-1/2" (14x89mm)	12/8.3	4518	Delrin	30	SW41Ti	TH-641
		4530	Kynar			
9/16x3-3/4" (14x95mm)	12/8.3	4536	Delrin	30	SW40Ti	---
		4537	Kynar			
1x3-1/2" (25x89mm)	38.5/25	4521	Delrin	23	SW28 SW32	AH629/36 Surespin 630/36
		4531	Kynar			

©2000-2015 Seton Scientific

- (1) Delrin is a registered trademark of E.I. DuPont de Nemours & Co. Delrin tubes must not be autoclaved.
- (2) Kynar is a registered trademark of Elf Atochem North America, Inc. Kynar tubes must not be autoclaved.
- (3) Tubes must only be used filled to within 2-3mm (1/8") of top (no partial fill) and used at 4 deg. C.

Continued on page 2

707-782-0900
707-782-0902 fax
sales@setonscientific.com
P.O. Box 750003
Petaluma, CA 94975 USA

OSR POLYPROPYLENE BOTTLES

Bottle Size (material)	Nominal Volume ml	Bottle Assy PN (1)	Bottle Only PN	Cap only PN	Plug Only PN (Delrin)	O-ring Only PN (Viton)	Maximum Speed KRPM (2)	For Beckman Rotor	For Thermo-Sorvall Rotor
5/8x3" (5a) 16x76mm (polypropylene)	10.4	4479 (alum cap)	5072	4418 (alum)	4500	4493	50	See (3)	T-1270,A-1256 T-865.1,T-875
5/8x3" (5a) 16x76mm (polypropylene)	10.4	4489 (Delrin cap)	5072	4499 (Delrin)	4500	4493	50	See (3)	T-1270,A-1256 T-865.1,T-875
5/8x3" (5b) 16x80mm (polypropylene)	9.1	5104 (alum cap)	4426 4425(7)	4424 (alum)	5081	4493	50	---	T-1270,A-1256 T-865.1,T-875
5/8x3" (5b) 16x80mm (polypropylene)	9.1	5105 (Delrin cap)	4426 4425(7)	5080 (Delrin)	5081	4493	50	---	T-1270,A-1256 T-865.1,T-875
1x3-1/2" (6a) 25x89mm (polypropylene)	25	5073 (alum cap)	5074	4419 (alum)	4500	5077	45	See (4)	---
1x3-1/2" (6a) 25x89mm (polypropylene)	25	5102 (Delrin cap)	5074	5103 (Delrin)	4500	5077	30	See (4)	---
1x3-1/2" (6b) 25x89mm (polypropylene)	23	5101 (alum cap)	5075	4501 (alum)	4500	5077	45	----	T-865 T-1250
1x3-1/2" (6b) 25x89mm (polypropylene)	23	5100 (Delrin cap)	5075	5076 (Delrin)	4500	5077	30	----	T-865 T-1250

©2000-2015 Seton

Scientific

(1) Assembly consists of one bottle, cap, plug and o-ring.

(2) Bottle must be used at full volume only regardless of the speed.

(3) Used in the following rotors: 90Ti, 80Ti, 75Ti, 70.1Ti, 65, 50Ti, 50 & 40.

(4) May be used in the following rotors: Types 70Ti, 60Ti, 55.2Ti, 50.2Ti, 45Ti, 42.1.

(5a) Beckman-type bottle. (5b) Sorvall-type bottle. Bottle/cap threads are different, parts are not interchangeable, but Beckman-type complete bottle assembly can be used in Sorvall rotors.

(6a) Beckman-type bottle. (6b) Sorvall-type bottle. Bottle/cap threads are the same but bottle lengths and cap/hole diameters are different. Parts are not interchangeable. Beckman-type to be used in Beckman rotors and Sorvall-type to be used in Sorvall rotors.

(7) PN 4425 is a box of 25 of PN 4426.