

SETON

EASY-SEAL TUBES Polyallomer & Polyclear Thin Wall *Ultra*

707-782-0900
707-782-0902 fax
sales@setonscientific.com
P.O. Box 750003
Petaluma, CA 94975 USA

EASY-SEAL POLYALLOMER & POLYCLEAR TUBES

Rev.3

Tube Size	Seton PN Box of 50	Tube Volume ml	Material	Use Crown Assy	For Beckman Rotor	For Thermo/Sorvall Rotor	Thermo Sorvall Tube PN	Thermo Sorvall Crown PN
7/16x1" (11x29mm)	6011	1.5	polyallomer	4264	TLA120.2	S150-AT S140-AT,S120-AT2	75000502	75000504
	8011		Polyclear				75000503	
1/2x1" (13x29mm)	6020	1.5	polyallomer	4267	TLA-100.3	S100-AT4 S100-AT6	75000512	75000511
	8020		Polyclear				---	
1/2x1-1/2" (13x44mm)	6021	3.5	polyallomer	4267	TLA-100.3	S100-AT4 S100-AT6	75000509	75000511
	8021		Polyclear				75000510	
1/2x2" (13x51mm)	6022	4.4	polyallomer	4267	TLA-100.4	S100-AT6	75000508	75000511
	8022		Polyclear				---	
1/2x2-1/2" (13x64mm)	6023	5.6	polyallomer	4171	50.3Ti,40.3 40.2,50.4Ti	TFT45.6	---	---
	8023		Polyclear				---	
5/8x2-7/8" (16x73mm)	6041	10	polyallomer	4017,4018,4019 4020,4063,4064	see (1)	T-1270,T-890, A-1256,TFT80.13	---	---
	8041		Polyclear				---	
1x3-1/8" (25x79mm)	6051	28.5	polyallomer Polyclear	4111	SW32, SW28	Surespin 630/36 AH-629/36	---	---
1x3-1/2" (25x89mm)	6052	33	polyallomer	4015,4016	see (2)	---	---	---
	8052		Polyclear				---	
1-1/2x4" (38x102mm)	6060(3)	82	Polyallomer	4202	45Ti	---	---	---
	8060(3)		Polyclear				---	

©2000-2016 Seton

Scientific

(1) Used in the following rotors: 90Ti, 80Ti, 75Ti, 70.1Ti, 65, 50Ti, Ty50 & Ty40.

(2) Used in the following rotors: 70Ti, 60Ti, 55.2Ti, 50.2Ti, and 42.1.

(3) PN 6060 and 8060 are boxes of 25 tubes.