

# SETON

## MICRO-ULTRA TUBES Thin Wall, Easy-Seal, Re-Seal, Thick Wall, Cone-Top, Ultra-Cone

707-782-0900  
707-782-0902 fax  
[sales@setonscientific.com](mailto:sales@setonscientific.com)  
P.O. Box 750003  
Petaluma, CA 94975 USA

MICRO-ULTRA TUBES-see the main page of each tube type for more details

Rev.4

Tube Size	Seton PN Box of 50	Tube Type (1)	Volume MI (2) Nom-Act	Crown/Adapter PN Seton(Beckman)	Crown/Adapter PN Seton(Sorvall)	For Beckman Rotor	For Thermo-Sorvall Rotor	Like Beckman Tube PN	Like Sorvall Tube PN
7/16x1" (11x25mm)	9011	Re-Seal CTPA	1.5-1.5	4246 (344636)	---	TLS55, TLA-120.2 MLA-130	S55-S, S150-AT S140-AT, S120-AT2	344624	75000501
	9013(3)			---	4255 (45306)	---	S150-AT, S140-AT S120-AT2	---	---
7/16x1" (11x29mm)	6011	ESPA	1.5-1.5	4264 (---)	4264 (75000504)	TLA120.2	S150-AT S140-AT, S120-AT2	---	75000502
	8011	ESPC		4264 (---)	4264 (75000504)			---	75000503
7/16x1-1/4" (11x32mm)	9012	Re-Seal CTPA	2.0-1.9	4246 (344674)	---	TLS55	S55-S	344625	75000459
	9014(3)			---	4255 (45306)	---	S150-AT, S140-AT S120-AT2	---	45246
7/16x1-1/4" (11x32mm)	5061(4)	UCPA(5)	1.4-1.4	4226 (---)	4226 (75000470)	TLS55	S55S	---	75000517
	7061(4)	UCPC(6)		---	---			75000516	
7/16x1-3/8" (11x34mm)	5011	OTPA	2.2-2.2	---	---	TLS55	S55S	347357	45240
	7011	OTPC		---	---			347356	45315
7/16x1-3/8" (11x34mm)	5007	TWPA	1.0-1.0	---	---	TLS55, TLA-120.2 MLA-130	S55-S, S150-AT S140-AT, S120-AT2	---	45998
	2007	TWPCB	2.5-2.5					---	45237
	7070	UltraThin	---					---	75000515
1/2x1" (13x25mm)	9024	Re-Seal CTPA	2.0-1.9	4247(355936)	4247(75000457)	MLS-50, TLA-100.3	S52-ST, S100-AT4	345829	75000458
	9026(3)			---	4256(45318)	---	S100-AT4	---	75000513
1/2x1" (13x29mm)	6020	ESPA	1.5-1.7	4267 (---)	4267 (75000511)	TLA-100.3	S100-AT4 S100-AT6	---	75000512
	8020	ESPC		---	---			---	
1/2x1-1/4" (13x32mm)	9025	Re-Seal CTPA	3.5-3.1	4247(355936)	---	MLS-50	S52-ST	349621	75000456
	9027(3)			---	4256(45318)	---	S100-AT4	---	45247
1/2x1-1/2" (13x44mm)	6021	ESPA	3.5-3.4	4267 (---)	4267 (75000511)	TLA-100.3	S100-AT4 S100-AT6	---	75000509
	8021	ESPC		---	---			---	75000510
1/2x1-1/2" (13x45mm)	9021	RSPA	3.9	(342883)	---	TLN-100	---	358980	---
1/2x2" (13x51mm)	9028(3)	RSPA	5.0-5.1	---	4257(45321)	---	S100-AT6	---	45248
1/2x2" (13x51mm)	9029	Re-Seal CTPA	5.0-5.1	4249 (362307)	---	TLA-100.4	---	362248	---
1/2x2" (13x51mm)	5024	TWPA	3.0-3.5	---	---	TLA-100.3 MLS-50	S100-AT4 S52-ST	349623	46102
	2001	TWPCB						349622	46100, 45238
1/2x2" (13x51mm)	5022	OTPA	5.0	---	---	MLS-50	S52-ST	326819	3127
	7022	OTPC	5.0					344057	46101
	7072	UltraThin	5.8					---	---
1/2x2" (13x51mm)	6022	ESPA	5.0-4.4	4267 (---)	4267 (75000511)	TLA-100.4	S100-AT6	---	75000508
	8022	ESPC		---	---			---	
1/2x2-1/4" (13x56mm)	5026	TWPA	3.2-3.2	---	---	TLA-100.4	---	362333	45598
	2010	TWPCB						362305	---
5/8x2" (16x53mm)	9042	RSPA	8.3	---	4491	---	S80-AT3	---	45597
	5047		5.3					---	45596
	2018		---					45595	
5/8x2-1/2" (16x64mm)	9040	Re-Seal CTPA	6.3-6.8	4248	---	MLA-80	---	345830	---

©2000-2016 Seton Scientific

- (1) *Abbreviations:* CTPA= Cone-Top Polyallomer; ESPA=Easy-Seal Polyallomer; ESPC=Easy-Seal Polyclear; Cone-Top Polyallomer; Re-Seal Polyallomer; TWPA=Thickwall Polyallomer; OTPA=Open-Top Polyallomer; OTPC=Open-Top Polyclear; TWPCB=Thickwall Polycarbonate; UCPA=UltraCone Polyallomer; UCPC=UltraCone Polyclear.
- (2) Nominal and actual volume for Seton tube.
- (3) Tubes have a shorter stem.
- (4) Requires UltraCone adapter PN 4226. To determine rMAX used in k-factor calculations, subtract the adapter thickness, 3.88mm.