

SETON

OPEN-TOP ULTRATHIN TUBES Very Thin Wall *Ultra*

707-782-0900
707-782-0902 fax
sales@setonscientific.com
P.O. Box 750003
Petaluma, CA 94975 USA

Thin wall tubes in swinging bucket rotors must be filled to 3mm of tube top!

ULTRATHIN TUBES

Rev. 2

Tube Size	Polyclear PN Box of 50	UltraCote PN Box of 50	Tube Capacity ml	Maximum Speed KRPM	For Beckman Rotor	For Thermo- Sorvall Rotor	For Kontron Rotor	Like Thermo- Sorvall Tube PN
7/16x1-3/8" (11x34mm)	7070	7070W	2.5	55	TLS55	RP55-S	---	75000515 (7070)
7/16x2-3/8" (11x60mm)	7071	7071W	4.7	60	SW60Ti SW56Ti	TH-660	TST60.4	
1/2x2" (13x51mm)	7072	7072W	5.8	65	SW65Ti, SW55Ti SW50.1, SW50, MLS-50	AH-650	TST55.5	
9/16x3-1/2" (14x89mm)	7073	7073W	12.8	41	SW41Ti	TH-641	---	
9/16x3-3/4" (14x95mm)	7074	7074W	15.7	41	SW40Ti	---	TST41.14	
5/8x2" (16x51mm)	7075	7075W	8.5	30	SW30.1	---	---	
5/8x4" (16x102mm)	7076	7076W	18.6	28	SW32.1 SW28.1Ti SW27.1Ti	AH-629/17 Surespin 630/17	TST28.17	
1x2" (25x51mm)	7077	7077W	15.1	30	SW30	AH-629/20	---	
1x3" (25x76mm)	7078	7078W	27.4	25	SW25.1	---	---	
1x3-1/2" (25x89mm)	7079	7079W	40	32	SW27, SW28 SW32	AH-629/36 T-1250 Surespin 630/36	TST28.38	

©2000-2015 Seton Scientific