

SETON

TUBES & BOTTLES LISTED BY ROTOR Beckman-Coulter

707-782-0900
707-782-0902 fax
sales@setonscientific.com
P.O. Box 750003
Petaluma, CA 94975 USA
rev.1

Only tube and bottle part numbers are listed here under each rotor currently manufactured by Beckman Coulter. With few exceptions the tubes and bottles require caps, crowns or spacers for sealing and support in the listed rotors. The corresponding caps, crowns and spacers are specified in the product page the part number link will direct you to.

* mfg. discontinued

* Only rotors presently manufactured (2020) are listed with tubes and bottles links. Beckman Coulter has discontinued the manufacture of many rotors including:

SW30, SW25.1, 60Ti, 75Ti, TY65, TY42.1, 55.2Ti, TY35, TY21, TY50, TY50.3, VTi80, NVT80, VTi65, 80Ty, SW65Ti, SW50.1, SW50.

If you are considering the use of a discontinued rotor please be sure the warranty has not expired and if a speed derating applies.

Please contact Seton Scientific Inc. if you have a rotor not listed and we can assist you with choosing a tube or bottle.

ABBREVIATIONS:

PCB= Polycarbonate
PA= Polyallomer
PP=Polypropylene

Fixed-Angle Rotors

90Ti, 70.1Ti, TY65*

5/8" (16 mm) Dia.

SSI PN	TYPE
2004	Thick Wall PCB
2015	PCB Bottle
3041	Re-Seal Polyclear
4281	Threaded Tube PCB
5041	Open-Top PA
5045	Thick Wall PA
5072	PP Bottle
6041	Easy-Seal PA
7041	Open-Top PC
8041	Easy-Seal PC
9040	Cone-Top PA
9041	Re-Seal PA

MLA-55, MLA-80

5/8" (16 mm) Dia.

SSI PN	TYPE
5044	Thick Wall PA
2003	Thick Wall PCB
9040	Cone-Top PA

70Ti, 50.2Ti, 60Ti*

1" (25 mm) Dia.

SSI PN	TYPE
2009	Thick Wall PCB
2016	PCB Bottle
3052	Re-Seal PC
4286	Threaded Tube PCB
5052	Open-Top PA
5055	Thick Wall PA
5074	PP Bottle
6052	Easy-Seal PA
7052	Open-Top PC
8052	Easy-Seal PC
9050	Cone-Top PA
9052	Re-Seal PA

Swinging Bucket Rotors

SW28Ti, SW32Ti

1" (25 mm) Dia.

SSI PN	TYPE
2009	Thick-Wall PCB
5052	Open-Top PA
5055	Thick Wall PA
5067	Ultra-Cone PA
6051	Easy-Seal PA
7052	Open-Top PC
7052W	Open-Top PC Wettable
7067	Ultra-Cone PC
7079	UltraThin
7079W	UltraThin Wettable
9050	Cone-Top PA
9051	Cone-Top PA

SW28.1Ti, SW32.1Ti

5/8" (16 mm) Dia.

SSI PN	TYPE
2005	Thick Wall PCB
5042	Open-Top PA
5046	Thick Wall PA
5066	Ultra-Cone PA
7042	Open-Top PC
7042W	Open-Top PC Wettable
7066	Ultra-Cone PC
7076	UltraThin
7076W	UltraThin Wettable
9040	Cone-Top PA
9056	Cone-Top PA

(Swinging Bucket Rotors continued, next page)

Vertical Rotors

VTi50

1" (25 mm) Dia.

SSI PN	TYPE
3052	Re-Seal PC
9052	Re-Seal PA

VTi65.1

5/8" (16 mm) Dia.

SSI PN	TYPE
3041	Re-Seal PC
9040	Cone-Top PA
9041	Re-Seal PA

VTi65.2, VTi90

1/2" (13 mm) Dia.

SSI PN	TYPE
3022	Re-Seal PC
9022	Re-Seal PA

Near Vertical Rotors

MLN80

5/8" (16 mm) Dia.

SSI PN	TYPE
9040	Cone-Top PA

NVT100, NVT90, NVT65.2

1/2" (13 mm) Dia.

SSI PN	TYPE
3022	Re-Seal PC
9022	

(Continued, next page)

Fixed-Angle Rotors (continued)45Ti
1-1/2" (38 mm) Dia.

SSI PN	TYPE
2011	Thick Wall PCB
4385	Threaded Tube PCB
5057	Thick Wall PA
5058	Thick Wall PA
5060	Open-Top PC
6060	Easy-Seal PA
7060	Open-Top PC
9060	Re-Seal PA

TLA-120.2, MLA-150, MLA-130
7/16" (11 mm) Dia.

SSI PN	TYPE
2007	Thick Wall PCB
5007	Thick Wall PA
6011	Easy-Seal PA
8011	Easy-Seal PC
9011	Cone-Top PA
9012	Cone-Top PA

100Ti, 50.4Ti, TLA-110, TLA-100.3
1/2" (13 mm) Dia.

SSI PN	TYPE
2001	Thick Wall PCB
5024	Thick Wall PA
9024	Cone-Top PA
9025	Cone-Top PA
9029	Re-Seal PA
9058	Re-Seal PA

Swinging-Bucket Rotors (continued)SW41
9/16" (14 mm) Dia.

SSI PN	TYPE
2013	Thick Wall PCB
5030	Open-Top PA
5032	Thick Wall PA
5064	Ultra-Cone PA
7030	Open-Top PC
7030W	Open-Top PC Wettable
7064	Ultra-Cone PC
7073	UltraThin
7073W	UltraThin Wettable
9032	Cone-Top PA
9055	Cone-Top PA

SW40Ti
9/16" (14 mm) Dia.

SSI PN	TYPE
2014	Thick Wall PCB
5031	Open-Top PA
5033	Thick Wall PA
5065	Ultra-Cone PA
7031	Open-Top PC
7031W	Open-Top PC Wettable
7065	Ultra-Cone PC
7074	UltraThin
7074W	UltraThin Wettable
9032	Cone-Top PA
9055	Cone-Top PA

SW55Ti, MLS-50
1/2" (13 mm) Dia.

SSI PN	TYPE
2001	Thick Wall PCB
5022	Open-Top PA
5024	Thick Wall PA
5063	Ultra-Cone PA
7022	Open-Top PC
7063	Ultra-Cone PC
7072	UltraThin
7072W	UltraThin Wettable
9024	Cone-Top PA
9025	Cone-Top PA
9058	Cone-Top PA

Near Vertical Rotors (continued)NVT65
1/2" (13 mm) Dia.

SSI PN	TYPE
3041	Re-Seal PC
9041	Re-Seal PA