









# SETON RE-SEAL CROWNS

## Ultra

707-782-0900  
 707-782-0902 fax  
[sales@setonscientific.com](mailto:sales@setonscientific.com)  
 P.O. Box 750003  
 Petaluma, CA 94975 USA

### ALUMINUM CROWNS

Rev.3

 <p><b>4065</b>                      +Aluminum,                      clear anodize                      +60 KRPM Max                      +Use Tubes:                      9053                      +For Rotors:                      70Ti,60Ti,55.2Ti,                      50.2Ti</p>	 <p><b>4068</b>                      +Aluminum,                      clear anodized                      +75 KRPM Max                      +Use Tubes:                      9041R                      +For Rotors:                      90Ti,80Ti,75Ti,                      70.1Ti</p>	 <p><b>4069</b>                      +Aluminum,                      Black anodized                      +65 KRPM Max                      +Use Tubes:                      9041R                      +For Rotors:                      Ty65,50Ti</p>
 <p><b>4491</b>                      +Aluminum,                      Black anodized                      +80 KRPM Max                      +Use Tubes:                      9042                      +For Rotors:                      S80-AT3,S58-A</p>	 <p><b>4257</b>                      +Aluminum,                      Black anodized                      +100 KRPM Max                      +Use Tubes:                      9028                      +For Rotors:                      S100-AT6</p>	 <p><b>4146</b>                      +Aluminum,                      clear anodize                      +60 KRPM Max                      +Use Tubes:                      9054                      +For Rotors:                      S55-A</p>
 <p><b>4492</b>                      +Aluminum,                      clear anodize                      +45 KRPM Max                      +Use Tubes:                      9060                      +For Rotors:                      45Ti</p>	 <p><b>4260</b>                      +Aluminum,                      Black anodized                      +58 KRPM Max                      +Use Tubes:                      9043                      +For Rotors:                      S58-A</p>	