

# SETON

## TUBES & BOTTLES LISTED BY ROTOR Thermo-Sorvall

707-782-0900  
707-782-0902 fax  
[sales@setonscientific.com](mailto:sales@setonscientific.com)  
P.O. Box 750003  
Petaluma, CA 94975 USA  
rev.1

**NOTES:**

Only tube and bottle part numbers are listed here under each rotor currently manufactured by Thermo-Sorvall. With few exceptions the tubes and bottles require caps, crowns or spacers for sealing and support in the listed rotors. The corresponding caps, crowns and spacers are specified in the product page the part number link will direct you to.

If you are considering the use of a discontinued rotor please be sure the warranty has not expired and if a speed de-rating applies.

Please contact Seton Scientific Inc. if you have a rotor not listed and we can assist you with choosing a tube or bottle.

**ABBREVIATIONS:**

PA= Polyallomer  
PC= Polyclear  
PCB= Polycarbonate  
PP= Polypropylene

### **Fixed-Angle Rotors**

T-8100, TFT-45.6  
1/2" (13 mm) Dia.

**SSI PN**  
2002  
5025  
9024  
9025  
9029  
9058

TFT 80.4  
7/16" (11 mm) Dia.

**SSI PN**  
2008  
4387  
5008  
5012  
9011  
9012  
9059

TFT 80.2  
5/16" (8 mm) Dia.

**SSI PN**  
5009

*(Fixed-Angle Rotors continued, next page)*

### **Swinging Bucket Rotors**

TH-660  
7/16" (11 mm) Dia.

**SSI PN**  
2008  
5008  
5010  
5062  
7010  
7010W  
7062  
7071  
7071W  
9011  
9012  
9059

AH-650  
1/2" (13 mm) Dia.

**SSI PN**  
2001  
5022  
5024  
5063  
7022  
7063  
7072  
7072W  
9024  
9025  
9058

*(Swinging Bucket Rotors continued, next page)*

### **Vertical Rotors**

StepSaver 70V6  
1/2" (13 mm) Dia.

**SSI PN**  
3022  
9022

StepSaver 65V13  
5/8" (16 mm) Dia.

**SSI PN**  
3041  
9041

### **Micro-Ultra Rotors**

**Fixed-Angle**  
S150-AT, S140-AT, S120-AT2  
7/16" (11 mm) Dia.

**SSI PN**  
2007  
5007  
6011  
8011  
9011  
9013  
9014

S110-AT, S100-AT4  
1/2" (13 mm) Dia.

**SSI PN**  
2010  
6020  
6021  
6022  
8021  
9028

*(Micro-Ultra Rotors continued, next page)*

**Fixed-Angle Rotors** (continued)

T-890, T-1270  
5/8" (16 mm) Dia.

**SSI PN**  
2004  
2015  
4425  
4426  
4494  
4495  
5045  
5072  
6041  
8041

T-865, T-1250  
1" (25 mm) Dia.

**SSI PN**  
5075  
5079

T-647.5  
1-1/2" (34 mm) Dia.

**SSI PN**  
5087  
5092

A-621  
61 mm Dia.

**SSI PN**  
4383

**Swinging Bucket Rotors** (continued)

TH-641  
9/16" (14 mm) Dia.

**SSI PN**  
2013  
5030  
5032  
5064  
7030  
7030W  
7064  
7073  
7073W  
9032  
9055

A-629/36, Surespin 630/36,  
Surespin 632/38.5  
1" (36 mm) Dia.

**SSI PN**  
2009  
5052  
5055  
5067  
6051  
7052  
7052W  
7067  
7079  
7079W  
9050  
9051

**Micro-Ultra Rotors, Fixed Angle** (continued)

S80-AT3  
16 mm Dia.

**SSI PN**  
2018  
5047  
9042

S58-A  
16 mm Dia.

**SSI PN**  
2021  
2022  
4426  
4495  
9043

S50-A  
25 mm Dia.

**SSI PN**  
2023  
2024  
9054

**Micro-Ultra Rotors  
Swinging Bucket**  
S55-S

7/16" (11 mm) Dia.

**SSI PN**  
2007  
5006  
5007  
5011  
5061  
7011  
7011W  
7061  
7070  
7070W  
9011  
9012

S52-ST  
1/2" (13 mm) Dia.

**SSI PN**  
2001  
5022  
5024  
5063  
7022  
7063  
7072  
7072W  
9024  
9025  
9058

(Micro-Ultra Rotors  
continued, next page)

**Micro-Ultra Rotors, Swinging Bucket**  
(continued)

S50-ST  
16 mm Dia.

**SSI PN**  
5048  
7048

**Micro-Ultra Rotors**  
**Vertical Rotors**

S120-VT  
7/16" (11 mm) Dia.

**SSI PN**  
9014