







# SETON

## CONE-TOP CROWNS Re-Seal *Ultra*

707-782-0900  
707-782-0902 fax  
[sales@setonscientific.com](mailto:sales@setonscientific.com)  
P.O. Box 750003  
Petaluma, CA 94975 USA

NORYL & DELRIN CROWNS

Rev.3

 <p><b>4216</b> +Noryl +70 KRPM Max +Use Tubes: 9050, 9051 +For Rotors: 70Ti,60Ti,55.2Ti, 50.2Ti,SW32,SW28</p>	 <p><b>4219</b> +Noryl +41 KRPM Max +Use Tubes: 9032, 9055 +For Rotors: SW41, SW40, TH-641</p>	 <p><b>4246</b> +Noryl +100 KRPM Max +Use Tubes: 9011, 9012 +For Rotors: TLS55, S55-S, TH660,TFT80.4</p>
 <p><b>4247</b> +Noryl +100 KRPM Max +Use Tubes: 9024, 9025, 9058 +For Rotors: SW55,100Ti,50.4Ti TFT45.6,AH-650</p>	 <p><b>4248</b> +Noryl +80 KRPM Max +Use Tubes: 9040, 9056 +For Rotors: 90Ti,80Ti,75Ti,65 50Ti,50,40,SW32.1 SW28.1</p>	 <p><b>4249</b> +White Delrin +100 KRPM Max +Use Tubes: 9029 +For Rotors: TLA-100.4</p>

©2000-2021 Seton Scientific